

MATERIAL SAFETY DATA SHEET

IDENTITY (AS USED ON LABEL AND LIST)

Safetec Antibiotic Ointment

Page 2 of 2

Date: September 21,2007

SECTION 4 - FIRE FIGHTING MEASURES

FLASH POINT (METHOD USED) Not available	FLAMMABLE LIMITS (% Volume in Air for Lowest Flashing Component) LEL: Not available UEL: Not available
--	---

EXTINGUISHING MEDIAUse extinguishing media appropriate for the surrounding fire. Use alcohol foam, dry chemical or CO₂.**SPECIAL FIRE FIGHTING PROCEDURES**

Wear full protective clothing.

UNUSUAL FIRE AND EXPLOSIVE HAZARDS

May emit toxic fumes.

SECTION 5 - PRECAUTIONS FOR SAFE HANDLING AND USE / LEAK PROCEDURES**STEPS TO BE TAKEN IF MATERIAL IS SPILLED OR RELEASED**

Assure slippery conditions are nullified. Add sand or dirt. Scoop up and remove. Wash area with soap and water. Dispose in accordance with federal, state, and local regulations.

WASTE DISPOSAL METHODS

Dispose of in accordance with Local, State, and Federal regulations. Products classified as non-hazardous may become hazardous waste upon contact with other products. Refer to "40 CFR Protection of Environment Parts 260-299" for complete waste disposal regulations. Consult your Local, State, or Federal Environmental Protection Agency before disposing of any chemicals.

PRECAUTIONS TO BE TAKEN IN HANDLING AND STORAGE

Keep containers closed when not in use. Do not handle or store near heat, sparks, flames or oxidizers.

SECTION 6 - EXPOSURE CONTROLS / PERSONAL PROTECTION

RESPIRATORY PROTECTION Not required	VENTILATION LOCAL EXHAUST: Recommended MECHANICAL (GENERAL): Not required
PROTECTIVE GLOVES Not required for normal use	EYE PROTECTION Not required for normal use
OTHER PROTECTIVE CLOTHING OR EQUIPMENT Eye wash stations should be nearby and ready for use.	WORK / HYGIENIC PRACTICES Practice safe work habits. Use according to label instructions.

SECTION 7 - PHYSICAL / CHEMICAL PROPERTIES

BOILING POINT 343°C (650 °F)	SPECIFIC GRAVITY (WATER=1) 0.855
VAPOR PRESSURE (mm Hg) No data	MELTING POINT 54°C (130 °F)
VAPOR DENSITY (AIR=1) No data	PH (1% solution in water) No data
SOLUBILITY IN WATER Insoluble	EVAPORATION RATE (IPA=1) No data
APPEARANCE AND ODOR White to off- white ointment	% VOLATILES (BY VOLUME) No data

SECTION 8 - STABILITY AND REACTIVITY

STABILITY UNSTABLE: STABLE: XXX	CONDITIONS TO AVOID Heat and oxidizing
INCOMPATIBILITIES (MATERIALS TO AVOID) Extreme heat, sparks, open flame and oxidizers.	
HAZARDOUS DECOMPOSITION OR BYPRODUCTS May give off fumes of carbon monoxide and carbon dioxide.	
HAZARDOUS POLYMERIZATION MAY OCCUR: WILL NOT OCCUR: XXX	CONDITIONS TO AVOID Temperatures above 300°C,

The information contained herein is believed to be accurate but is not warranted to be so. Users are advised to confirm in advance of need that information is current, applicable and suited to the circumstances of use. Vendor assumes no responsibility for injury to vendee or third persons proximately caused by the material if reasonable safety procedures are not adhered to as stipulated in the data sheet. Furthermore, vendor assumes no responsibility for injury caused by abnormal use of this material even if reasonable safety procedures are followed.